



Laerdal

helping save lives

Little Junior QCPR



Bystander CPR saves lives

The Little Junior QCPR manikin provides realistic child CPR training to support the community's role in improving the outcomes of cardiac arrests.

Little Junior QCPR makes it easy for instructors to focus on teaching quality CPR to every student, while adding a fun and engaging learning experience.





Fun and engaging

With objective feedback and scoring, learners are motivated not only to compete with one another, but also to try to beat their own best score. The friendly competition element makes training more engaging and interactive and has been demonstrated to improve motivation.

Efficient teaching for effective learning



Little Junior QCPR allows an instructor to monitor up to six learners at the same time from the palm of their hand. The QCPR Instructor app provides an easy to view performance summary that highlights which students may require further guidance. At a glance you can see who is not compressing deeply enough, who is not compressing at the correct rate and who is overventilating. Students can monitor their own performance with the QCPR Learner app or the tethered SkillGuide device.

QCPR measurement and feedback technology enables:

- Real-time feedback on compressions and ventilations
- Post-training debrief, including overall scores and tips for improvement
- Details on compression release, depth and rate, ventilation volume, and number of compressions, ventilations and cycles




QCPR Learner App
QCPR Instructor App



QCPR Learner App

All existing Little Junior manikins can be quickly and easily upgraded to QCPR with the Little Junior QCPR upgrade kit.



Laerdal
helping save lives